

BUSTER Extreme Dog Muzzle

- Easy to Fit
- Tough and durable, yet flexible and soft
- Variable collar adjustment
- Safe and strong dog muzzle with strategic addition of struts at front nose area to eliminate all possibilities of biting
- Dog friendliness in terms of feeding (nibble and water) and panting

The BUSTER muzzle is a comfortable dog muzzle designed to protect against bites.

The safety strapping ensures that the dog muzzle will always remain in place and it features two additional points of secure attachment. The loop at the bottom is designed to attach to your dog's regular collar.

The muzzle has a secure metal buckle with pre-holed webbing to ensure easy adjustment and a good fit.

The BUSTER Extreme Dog Muzzle is made from Injection moulded thermo plastic rubber, which is extremely tough and durable.

The neoprene padding inside of the webbing straps ensures that it is snug and comfortable for the dog to wear.

The dog can pant, drink and be given treats to reward a more placid attitude to life.



Tough and durable material for security purpose



Allows drinking, panting and vomiting



Inhibits biting, chewing and barking



Difficult for dog to dislodge




Sizing

The best way to choose the model size is to use the weight/breed guidelines, as long as your dog is not over or underweight.

If your dog is not at an ideal weight, or is between sizes, we recommend that you use these measuring instructions:

- Use a soft measuring tape, measure the length of your dog's snout from the end of his nose to the point where the snout starts to slope upwards to the eyes as shown on the drawing
- Then, measure around the widest part of your dog's snout while his mouth is closed



Weight/Breed Guideline:

Cat. No	Size	Approx. Weight of dog (kg/lbs)	Example Breeds (Avg. Sized adults)	Snout Length (cm / inch)	Snout Circumference (cm / inch)	Adjustable Neck Circumference (cm/ Inch)
272351	Size 1	4.5 to 7 kg 10 to 15 lbs	Bischoon Frise, Dachshund, Miniature Pinscher, Large Chihuahua, Jack Russell	5.5 cm / 2.1 in	21.5 cm / 8.4 in	20 cm / 7.8 in to 39 cm / 15.3 in
272352	Size 2	5.5 to 11 kg 12 to 25 lbs	Border Terrier, Jack Russell, Cocker Spaniel, Small Beagle, Miniature Poodle, Sheltie	6.5 cm / 2.5 in	27 cm / 10.6 in	21 cm / 8.2 in to 44 cm / 17.3 in
272353	Size 3	11 to 20 kg 25 to 45 lbs	Border Collie, Large Beagle, Corgi, Vizsla, Standard Poodle	7 cm / 2.7 in	30 cm / 11.8 in	27 cm / 10.6 in to 52 cm / 20.4 in
272354	Size 4	18 to 29 kg 40 to 65 lbs	Springer Spaniel, Dalmatian, Pointers, Small Pit Bull, Cattle Dog	8 cm / 3.1 in	34 cm / 13.3 in	29 cm / 11.4 in to 58 cm / 22.8 in
272355	Size 5	27 to 40 kg 60 to 90 lbs	Labrador, Golden Retriever, Doberman, German Shepherd, Greyhound, Weimaraner, Husky, Large Pit Bull, Boxer	9.2 cm / 3.6 in	39 cm / 15.3 in	34 cm / 13.3 in to 60 cm / 23.6 in
272356	Size 6	36 to 68 kg 80 to 150 lbs	Great Dane, Mastiff, Newfoundland, Rottweiler, Giant Schnauzer, X-Large Pit Bull, Large German Shepherd	10.5 cm / 4.1 in	43 cm / 16.9 in	36 cm / 14.1 in to 71 cm / 27.9 in

The BUSTER Extreme Dog Muzzle is washable and can be sterilized. Do not use if damaged. Do not leave your dog muzzled without supervision.



Soft neoprene padding for extra comfort



Overhead strap connects directly to the rear strap



Allows panting and treating



Made in China
For veterinary use only

