






# BOVIVET Reinforced Disposable Needles

We offer a range of high quality reinforced disposable injection needles.

- A range of sizes specially designed for the veterinary profession
- Thickened needle hub for extra strength and to avoid break of needle hub
- Luer lock
- Non-pyrogenic
- The needle hub strength has been tested at 100-110N
- Made of high-quality stainless steel and polypropylene
- Sterilized by ethylene oxide
- Paper-film peel pack
- 100/pk

Cat. No	Ø mm x L mm	Gauge	Colour
112455	1.25x38mm	18Gx1½"	
112457	1.5x50mm	17Gx2"	
112458	1.6x25mm	16Gx1"	
112459	1.6x38mm	16Gx1½"	
112460	2.1x38mm	14Gx1½"	
112461	2.1x60mm	14Gx2½"	
112462	2.1x80mm	14Gx3 1/4"	