

**Mainatvet**  
by **KRUUSE**

# Canine Garments



## KRUUSE Abdominal Dog Belt


Recovery bandage wraps are made in polymer foam and are easily fitted and adjusted with velcro bands. "Trihook" extra-strong Velcro fixing bands which doesn't loose strenght in the worse working conditions (dust, dirt, hair, etc.). Wounds are covered with monofilm pad 1 cm thick which forms an air cushion to protect the wound.

All KRUUSE Abdominal dog belt have been tested on clinical cases in veterinary hospitals.

- Easily fitted
- Adjustable closures
- Breathable fabric
- Machine washable
- Maximum protection for the immediate post-operative period

### KRUUSE Abdominal female dog belt

Cat. No	Description	Size	Distance between legs
165400	KRUUSE Abdominal female dog belt, small breed	XS	12 – 14 cm
165401	KRUUSE Abdominal female dog belt, small breed	S	15 – 17 cm
165402	KRUUSE Abdominal female dog belt, small breed	M	18 – 20 cm
165403	KRUUSE Abdominal female dog belt, small breed	L	21 – 23 cm
165404	KRUUSE Abdominal female dog belt, medium breed	S	24 – 27 cm
165405	KRUUSE Abdominal female dog belt, medium breed	M	28 – 34 cm
165406	KRUUSE Abdominal female dog belt, medium breed	L	35 – 39 cm
165407	KRUUSE Abdominal female dog belt, large breed	S	40 – 45 cm
165408	KRUUSE Abdominal female dog belt, large breed	M	46 – 52 cm



How to measure the size

### KRUUSE Abdominal male dog belt

Cat. No	Description	Size	Distance between legs
165409	KRUUSE Abdominal male dog belt, small breed	S	15 – 17 cm
165410	KRUUSE Abdominal male dog belt, small breed	M	18 – 20 cm
165411	KRUUSE Abdominal male dog belt, small breed	L	21 – 23 cm
165412	KRUUSE Abdominal male dog belt, medium breed	S	24 – 27 cm
165413	KRUUSE Abdominal male dog belt, medium breed	M	28 – 34 cm
165414	KRUUSE Abdominal male dog belt, medium breed	L	35 – 39 cm
165415	KRUUSE Abdominal male dog belt, large breed	S	40 – 45 cm
165416	KRUUSE Abdominal male dog belt, large breed	M	46 – 52 cm



Scan QR code to see video

