

BUSTER Premium Dog Collar

The greatest innovation in e-collars since the BUSTER Collar!

*The collar is used as protection after surgery
The collar prevents the dog from reaching surgical wounds by licking, biting or chewing*

- **OPTIMAL COMPLIANCE**
- **QUICK APPLICATION**
- **SNUG FIT**

The BUSTER Premium Dog Collar is perfect for:

- Post operative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin lesions
- Skin disease



Scan the QR Code to watch the video

UNIQUE SELLING POINTS

From minutes to seconds

The time of assembling and application of the collar has been reduced from minutes to seconds

Great comfort

The banding around the outer edge will not mark or scrape furniture and it serves as a sound-absorbing feature to the dog

Safety and perfect fit

The collar is provided with extra neoprene and Velcro closing for safety and perfect fit

Soft, flexible and comfortable perforated neoprene

The collar is provided with soft, flexible and comfortable perforated neoprene

Avoid skin irritation

The banding around the inner edge provides better comfort and prevents skin irritation

Assembled collar ready to use

The collar is easy to fasten. You simply insert the dog's head into the collar and close the Velcro

Easy handling

The collar has a safe and easy Velcro closing

No need for fastening gauze strap or collar

The close-fitting perforated neoprene make the use of collar or fastening gauze straps unnecessary in the cage after surgery

Anatomic cut

The anatomic cut of the PP material follows the curve of the dog's neck and provides comfort to the dog during sleep or rest. The cut also eases the application for the care taker

Thought-through design ready to use



Guideline for Sizes

The BUSTER Premium Dog Collar is available in 5 sizes.

Cat. No	Packaging	Size	Body weight	Suitable for the following breeds	Dimensions of box:	Dimensions of master carton:
273528	5/pk	XX-small	Breeds from 1-3 kg	Chihuahua, Japanese Chin, Maltese, Pomeranian, Japanese Terrier, Yorkshire Terrier	D = 19 cm W = 19 cm H = 22.7 cm	D = 60 cm W = 60 cm H = 66 cm
273529	5/pk	X-small	Breeds from 3-6 kg	Bichon Frisé, Coton de Tulear, Papillon Dogs, American Hairless Terrier, Shih Tzu, Tibetan Spaniel, English Toy Spaniel	D = 21 cm W = 21 cm H = 23.2 cm	D = 75 cm W = 75 cm H = 75 cm
273531	5/pk	Small	Breeds 6-10 kg	Cairn Terrier, Miniature Schnauzer, Jack Russel Terrier, Dachshund, Poodle, Border Terrier, Cavalier King Charles Spaniel	D = 30 cm W = 30 cm H = 24 cm	D = 62 cm W = 62 cm H = 50 cm
273532	5/pk	Medium	Breeds 10-20 kg	Border Collie, Finnish Lapphund, Beagle, Staffordshire Bull Terrier, English Springer Spaniel, Cocker Spaniel	D = 39 cm W = 39 cm H = 29 cm	D = 80 cm W = 40 cm H = 60 cm
273533	5/pk	Large	Breeds 20-35 kg	German Shepherd, Labrador Retriever, Golden Retriever, Wirehaired Pointer, Boxer, Flatcoated Retriever	D = 39 cm W = 39 cm H = 29 cm	D = 50 cm W = 50 cm H = 71 cm
273534	5/pk		Starter Pack: 1 x XX-Small 1 x X-Small 1 x Small 1 x Medium 1 x Large		D = 47 cm W = 47 cm H = 31 cm	D = 48 cm W = 48 cm H = 65 cm

- **BETTER COMFORT**
- **SAFE AND PERFECT FIT**
- **EASY TO FASTEN**
- **READY TO USE**



*Easy stressfree handling
Great comfort*

Instructions for Use



Step 1
■ Insert the dog's head into collar



Step 2
■ Close the Velcro strap